





studiostand
junior model



STUDIO STAND Junior model

The ideal stand for those who have a small studio. Moreover, it is very suitable for use as a supplementary stand. The quality is identical to that of the other Cambo stands.

A geared centre column is provided. This centre column can be racked upwards, and an extra extension tube can be pulled out by hand to obtain an increase in height. At the front of the base the stand is provided with two castors and the back with a rubber buffer. The stand can be moved simply by lifting it slightly by means of handle G, thus the rubber buffer is lifted from the ground and the stand can be moved resting on the two front castors. The Junior studio stand can be supplied in two lengths, viz of 80 and 100 cm.

Order code for 80 cm stand: BST I
Order code for 100 cm stand: BST II

- A. Cast iron base
- B. Tube of 12 cm ϕ
- C. Minimum height of BST I 80 cm and BST II 100 cm
- D. Maximum height with racked up centre column and extended centre column 185 cm for BST I and 205 cm for BST II
- E. Aluminium geared centre column with extension tube
- F. Accessory tray size 30 x 40 cm
- G. Stand lifting handle
- H. Geared centre column winding handle

Weight BST I and BST II : 38 kgs

EXCLUSIVE U.S. DISTRIBUTOR

BROOKS®

BURLEIGH BROOKS OPTICS, INC.

44 BURLEWS COURT HACKENSACK N.J. 07601